

# Compreflex Reduce THIGH

Patient Name:	Contact Name:
Account Name:	Account #:
Account Phone #:	P.O. #:
Ship Name & Address:	Date:

## Product Information

Product includes one Compreflex Reduce Thigh wrap, one pair of Cotton Liners, and one Belt and Hip attachment.

<input type="checkbox"/> Right Thigh	Size:	Item #:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:
<input type="checkbox"/> Left Thigh	Size:	Item #:
Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige	Length:	Quantity:

## Companion Products (sold separately)

<input type="checkbox"/> Compreflex Reduce Knee	Color: <input type="checkbox"/> Black <input type="checkbox"/> Beige
Size:	Item #:      Quantity:
Size:	Item #:      Quantity:

## Measuring Instructions

### Compreflex Reduce Thigh:

Measure circumferences where indicated on the measurement chart, beginning 5cm above the top of patella and continue 23cm up the thigh, just below the gluteal fold. Record in the corresponding lines D and F. Measure length from the lateral aspect top of patella to the gluteal fold, and record in box G1.

### Compreflex Reduce Knee:

Measure circumferences where indicated, beginning 5cm below the lateral aspect fibular head. Continue 5cm above the lateral aspect top of patella. Record in the corresponding lines C1 and D.

## Circumference

Left      Right

F \_\_\_\_\_


D \_\_\_\_\_

C1 \_\_\_\_\_

**G1 Length**  
 Regular: 20–24.9cm  
 Tall: 25–30cm

*Top of patella to gluteal fold*

G1 \_\_\_\_\_



## Compreflex Reduce Thigh Size Chart

	Range 1	Range 3
F	<80cm	<125cm
D	<65cm	<110cm

		Black	
Left	Regular	2401-TCR-L	2403-TCR-L
	Tall	2401-TCT-L	2403-TCT-L
Right	Regular	2401-TCR-R	2403-TCR-R
	Tall	2401-TCT-R	2403-TCT-R
		Beige	
Left	Regular	2411-TCR-L	2413-TCR-L
	Tall	2411-TCT-L	2413-TCT-L
Right	Regular	2411-TCR-R	2413-TCR-R
	Tall	2411-TCT-R	2413-TCT-R

## Compreflex Reduce Knee Size Chart

	Range 1	Range 3
D	<75cm	<115cm
C1	<65cm	<90cm

		Black	
		2401-KP	2403-KP
		Beige	
		2411-KP	2413-KP

